



EFFECTIVE AGAINST
Staphylococcus Aureus (Staph),
Methicillin-resistant
Staphylococcus Aureus (MRSA),
and Vancomycin Resistant
Enterococcus Faecalis (VRE)

**THIS PATENTED PRODUCT
KILLS NOROVIRUS**
including murine norovirus
and Norwalk-like virus
(Feline calicivirus) in 30 seconds

This product is approved for use on Medical Equipment Surfaces. See inside label for details.

DIRECTIONS FOR USE: It is a violation of Federal Law to use this product in a manner inconsistent with its labeling. *This product is not to be used as a terminal sterilant/high level disinfectant on any surface or instrument that (1) is introduced directly into the human body, either into or in contact with the bloodstream or normally sterile areas of the body, or (2) contacts intact mucous membranes but which does not ordinarily penetrate the blood barrier or otherwise enter normally sterile areas of the body. This product may be used to pre-clean or decontaminate critical or semi-critical medical devices prior to sterilization or high-level disinfection.*

DISINFECTION DIRECTIONS
TO DISINFECT HARD, NONPOROUS SURFACES: Spray this product 6-8 inches from surface until thoroughly wet. Treated surfaces must remain wet for appropriate contact time listed on this label. Wipe or allow to air dry. No rinse required. *For heavily soiled surfaces, precleaning is required.*

TO DISINFECT NON CRITICAL PRECLEANED INSTRUMENTS: Instruments must be thoroughly cleaned to remove excess organic debris, rinsed and dried. Spray all surfaces of instruments with this product until thoroughly wet. Let stand for 1 minute. To kill TB let stand for 5 minutes. Wipe dry with a clean damp cloth or paper towel.

SANITIZATION DIRECTIONS
TO SPOT SANITIZE SOFT SURFACES*: Spray this product 6-8 inches from soft surface until wet. DO NOT SATURATE. Let stand for 20 seconds. Allow to air dry. *Heavy soils must be removed prior to sanitizing.*

SPECIAL INSTRUCTIONS FOR CLEANING AND DECONTAMINATION AGAINST HIV-1, HBV and HCV ON SURFACES/OBJECTS SOILED WITH BLOOD OR BODY FLUIDS: This product kills HIV, HBV AND HCV on precleaned environmental surfaces/objects previously soiled with blood/body fluids in healthcare settings or other settings in which there is an expected likelihood of soiling of inanimate surfaces/objects with blood/body fluids and in which the surfaces/objects likely to be soiled with blood/body fluids can be associated with the potential for transmission of Human Immunodeficiency Virus Type 1 (HIV 1) (associated with AIDS), Hepatitis B Virus (HBV), and Hepatitis C Virus (HCV)

Personal Protection: When handling items soiled with blood or body fluids use disposable impervious gloves, gowns, masks and eye coverings.

Cleaning Procedures: Blood and body fluids must be thoroughly cleaned from surfaces and objects before application of this product.

Contact Time: Spray on surface until completely wet and let stand 30 seconds. Note contact time for other microorganisms on this label may vary.

Infectious Materials Disposal: Blood and other body fluids should be autoclaved and disposed of according to local regulations for infectious disposal.

EPA Reg. No. 84368-1-84150
EPA Est. No. 9019-OH-001 (MA), 39754-WI-001 (CR)
Beginning of batch code indicates Est. No.
Patent Nos. 7,754,766; 8,110,604; 8,143,309

Distributed by: GOJO Industries, Inc. P.O. Box 991, Akron, OH 44309
Ingredient Information/Questions? Tel: 1-888-4-PURELL
©2016. GOJO Industries, Inc., All rights reserved. Made in U.S.A.

Reorder No. 3340 3340-641-B



TO REFILL: Remove trigger from empty bottle. Unscrew cap on refill and pour contents directly into empty bottle. Replace trigger sprayer and use as you normally would.

STORAGE AND DISPOSAL:

Store this product in a cool, dry area. Do not reuse or refill this container except as described in the directions for use. Refill only with this product. Offer for recycling or place in trash.

Physical or Chemical Hazards
Combustible. Do not use or store near heat or open flame.

BACTERIA:

30 second contact time

- Campylobacter jejuni*
- Escherichia coli (E. coli)*
- Klebsiella pneumoniae*
- Listeria monocytogenes (Listeria)*
- [†]Methicillin-resistant *Staphylococcus aureus (MRSA)*
- Salmonella enterica (Salmonella)*
- Streptococcus pneumoniae (Strep)*
- Streptococcus pyogenes (Strep)*
- [†]Vancomycin Resistant *Enterococcus faecalis (VRE)*
- Vibrio vulnificus*
- Yersinia enterocolitica*

1 minute contact time

- Pseudomonas aeruginosa*
- Staphylococcus aureus (Staph)*

Mold, Mildew & Fungi

- Aspergillus niger (Mold)* 5 minutes
- Candida albicans* 30 seconds
- Trichophyton mentagrophytes* 30 seconds

Mycobacterium

- Mycobacterium bovis var. BCG (TB)* 5 minutes
- [†]**Viruses Enveloped**
- [†]2009-H1N1 Influenza A Virus (H1N1) 30 seconds
- Influenza A virus (Flu Virus) 30 seconds

[†]**Viruses Non-Enveloped**

- Hepatitis A virus (HAV) 1 minute
- [†]Murine norovirus (Norovirus) 30 seconds
- [†]Feline Calicivirus (as surrogate for Norovirus) 30 seconds
- Polio Type 1 virus 30 seconds
- Rhinovirus (a cause of the common cold) 30 seconds

[†]**Bloodborne Pathogens**

- [†]Human Hepatitis B virus (HBV) 30 seconds
- [†]Human Hepatitis C virus (HCV) 30 seconds
- [†]Human immunodeficiency virus Type 1 (HIV-1) 30 seconds
- Soft Surface Sanitization Pathogens**
- Klebsiella pneumoniae* 20 seconds
- Staphylococcus aureus (Staph)* 20 seconds

Use Sites:

- Healthcare facilities
- Physician offices
- Dental offices
- Nursing homes
- Hospitals
- Nurseries
- Neonatal units
- Emergency rooms
- Ambulance
- Blood banks
- Dialysis centers
- Hospices
- Labs
- Medical research facilities
- Patient care rooms
- Exam rooms
- Physical therapy rooms
- Radiology rooms
- X-Ray labs
- Anesthesia rooms
- Surgery rooms
- Operating rooms

Medical Equipment Surfaces:

- External ultrasound transducers and probes
- MRI machines
- Chemotherapy hoods
- CAT Scans
- Operating tables
- Examination tables
- Dental chairs
- Stretchers

High Touch Surfaces:

- Bed frames & railings
- IV Poles
- Wheelchairs
- Nurse call buttons
- Cabinet handles
- Keyboards
- Light switches
- Doorknobs
- Tables
- Chairs

Hard Non-Porous Surface Materials:

- Acrylic
- Chrome
- Glazed ceramic tile
- Fiberglass
- Finished wood**
- Glass
- Laminated surfaces
- Metal
- Plastic
- Plexiglass®
- Porcelain
- Stainless steel
- Vinyl

***Soft Surfaces:**

- Hospital gowns
- Privacy drapes
- Scrub suits
- Lab coats
- Towels
- Linens
- Mattresses
- Couches
- Chairs
- Upholstered furniture

**Not recommended for use on alcohol soluble wood finishes. Test in small inconspicuous area before broad application.